

PETAPLUS®

TEETH MOULD CHART

ACRYLIC



COMPOSITE



SQS
ISO 13485

C€1250

A7
PARTNER AG

SWITZERLAND

PETAPLUS® L

Acrylic Teeth

●

11

1 39,8
2 8,7
3 7,1

■

14

1 40,3
2 8,0
3 7,4

●

17

1 42,3
2 8,0
3 7,7

●

20 A

1 43,0
2 9,0
3 7,9

▲

23 A

1 43,1
2 9,8
3 7,8

▲

26

1 43,5
2 9,4
3 7,9

■

29

1 44,1
2 9,3
3 8,1

●

32

1 44,5
2 9,5
3 8,2

■

35

1 45,3
2 9,5
3 8,2

▲

38 A

1 45,3
2 9,6
3 8,3

●

41

1 43,6
2 10,3
3 8,2

▲

44

1 46,5
2 10,7
3 8,7

■

47

1 49,3
2 12,0
3 9,5

■

51

1 33,3
2 7,6
3 5,2

■

52

1 32,0
2 7,9
3 4,8

■

53

1 32,3
2 9,0
3 4,9

■

54

1 32,7
2 9,0
3 4,8

■

55 A

1 34,8
2 8,7
3 4,8

■

56 A

1 35,8
2 8,6
3 5,4

■

57

1 35,5
2 8,7
3 5,5

■

58

1 34,7
2 9,0
3 5,1

■

59

1 36,0
2 10,3
3 5,5

■

60

1 37,2
2 11,0
3 5,6

62

1: 29,0 2: 32,0 3: 8,2 4: 8,9

65

1: 29,8 2: 31,2 3: 8,8 4: 7,8

67 A

1: 29,9 2: 33,0 3: 8,3 4: 8,9

68 A


1: 29,7 2: 33,0 3: 8,7 4: 8,7

69

1: 31,8 2: 35,6 3: 8,6 4: 9,3

78

1: 31,4 2: 33,3 3: 9,0 4: 9,0





COMBINATION CHART

11	51	62
14	52	62
17	53	62
20A	54	65
23A	55A	67A - 68A
26	54	65
29	57	65 - 67A
32	58	68A - 69
35	58	68A - 69
38A	56A	68A - 69
41	59	69 - 78
44	60	78
47	60	78





G1	22	29
G2	25	31
G3	25-34	31-33
G4	37	33-35
G5	25-34	33
G6	37	35
G7	37	35-39
F1	38	35-39
F2	38	35-39
F3	38	35-39
F4	67	39
F5	74	39-41
F6	88	41
F7	58	41
H1	25	31
H2	34	33
H3	44	35-39
H4	44	35-39



PETAPLUS[®] H

Acrylic Teeth



K1

1 40,4
2 8,3
3 7,4

K3

1 40,0
2 8,8
3 7,4

K4

1 41,6
2 9,6
3 7,9

K5

1 42,7
2 8,4
3 8,3

K6

1 41,9
2 9,5
3 7,8

M1

1 43,7
2 9,6
3 8,1

M2

1 43,5
2 8,9
3 8,1

M3

1 44,0
2 10,0
3 8,3

M4 A

1 45,0
2 9,7
3 8,3

M5

1 43,6
2 10,0
3 8,2

L1

1 47,1
2 9,2
3 9,0

L2

1 47,3
2 9,3
3 8,9

L3

1 46,4
2 10,5
3 8,6

L4

1 47,3
2 10,5
3 8,8

L5 A

1 49,0
2 12,4
3 9,3

E1

1 31,7
2 7,6
3 4,5

E2

1 33,7
2 9,3
3 5,0

E3

1 34,5
2 9,3
3 5,4

E4

1 33,7
2 7,6
3 4,9

E5

1 37,0
2 9,2
3 5,5

E6

1 39,7
2 9,1
3 6,3

E7 A

1 38,2
2 12,0
3 5,5

N1

1: 29,2 2: 29,9 3: 8,2 4: 7,7

N2

1: 29,0 2: 31,0 3: 8,5 4: 8,0

N3

1: 30,0 2: 33,0 3: 8,0 4: 8,4

N4


1: 31,2 2: 33,0 3: 8,2 4: 8,6

N5

1: 30,5 2: 33,4 3: 8,9 4: 8,7

N6

1: 34,2 2: 35,5 3: 9,4 4: 9,4



COMBINATION CHART



K1	E1	N1
K3	E1	N1 - N2
K4	E2	N2
K5	E2 - E4	N2 - N3
K6	E2	N2 - N3
M1	E2	N2 - N3
M2	E2 - E4	N2 - N3
M3	E3 - E4	N3 - N4
M4A	E3	N4 - N5
M5	E3 - E4	N3 - N4
L1	E3 - E6	N4 - N5 - N6
L2	E3 - E6	N5 - N6
L3	E5	N5 - N6
L4	E5	N5 - N6
L5A	E7A	N6



PETAPLUS[®] C

Composite Teeth



K1

1 40,4
2 8,3
3 7,4

K3

1 40,0
2 8,8
3 7,4

K4

1 41,6
2 9,6
3 7,9

K5

1 42,7
2 8,4
3 8,3

K6

1 41,9
2 9,5
3 7,8

M1

1 43,7
2 9,6
3 8,1

M2

1 43,5
2 8,9
3 8,1

M3

1 44,0
2 10,0
3 8,3

M4
A

1 45,0
2 9,7
3 8,3

M5

1 43,6
2 10,0
3 8,2

L1

1 47,1
2 9,2
3 9,0

L2

1 47,3
2 9,3
3 8,9

L3

1 46,4
2 10,5
3 8,6

L4

1 47,3
2 10,5
3 8,8

L5
A

1 49,0
2 12,4
3 9,3

E1

1 31,7
2 7,6
3 4,5

E2

1 33,7
2 9,3
3 5,0

E3

1 34,5
2 9,3
3 5,4

E4

1 33,7
2 7,6
3 4,9

E5

1 37,0
2 9,2
3 5,5

E6

1 39,7
2 9,1
3 6,3

E7
A

1 38,2
2 12,0
3 5,5

N1

1: 29,2 2: 29,9 3: 8,2 4: 7,7

N2

1: 29,0 2: 31,0 3: 8,5 4: 8,0

N3

1: 30,0 2: 33,0 3: 8,0 4: 8,4

N4


1: 31,2 2: 33,0 3: 8,2 4: 8,6

N5

1: 30,5 2: 33,4 3: 8,9 4: 8,7

N6

1: 34,2 2: 35,5 3: 9,4 4: 9,4



COMBINATION CHART



K1	E1	N1
K3	E1	N1 - N2
K4	E2	N2
K5	E2 - E4	N2 - N3
K6	E2	N2 - N3
M1	E2	N2 - N3
M2	E2 - E4	N2 - N3
M3	E3 - E4	N3 - N4
M4A	E3	N4 - N5
M5	E3 - E4	N3 - N4
L1	E3 - E6	N4 - N5 - N6
L2	E3 - E6	N5 - N6
L3	E5	N5 - N6
L4	E5	N5 - N6
L5A	E7A	N6





EC Certificate

SQS as a conformity assessment body identification number 1250 herewith certifies the organisation

AZ & PARTNER AG

**Dorfchärn
6243 Egolzwil
Switzerland**

the use of a quality assurance system in its manufacturing and distribution which fulfills the requirements set out in:

ANNEX V

Directive 93/42/EEC

This approval is based on the result of the report dated April 1, 2020.

The scope of validity covers the products

**Dental alloys
Artificial teeth
Zircon blanks**

The following CE label can be applied to these products mentioned in the Appendix of this certificate

CE 1250

A condition for the validity of this certificate is a regular examination in accordance with Annex V.4 of the Directive 93/42/EEC.

Reg. no. 41686

Validity 29.04.2020 – 25.05.2024
Issue 29.04.2020

Approved Medical Responsible
29.04.2020

F. Müller, CEO SQS

D. Taddeo, Medical Responsible



Swiss Association for Quality
and Management Systems (SQS)
Bernstrasse 103, 3052 Zollikofen, Switzerland



ANNEX V

Appendix to the EC Certificate

Directive 93/42/EEC

This Appendix is valid only in connection with the following certificate:

Registration Number 41686

Validity from April 29, 2020 up to and including May 25, 2024

This approval includes the following Medical Device/s:

Class IIa
METAPLUS dental alloys
PETAPLUS artificial teeth
ZIRKOPLUS zircon blanks

Appendix Issue: April 29, 2020

Reg. no. 41686

Validity 29.04.2020 – 25.05.2024
Issue 29.04.2020

Approved Medical Responsible
29.04.2020

F. Müller, CEO SQS

D. Taddeo, Medical Responsible



Swiss Association for Quality
and Management Systems (SQS)
Bernstrasse 103, 3052 Zollikofen, Switzerland





PARTNER AG

SWITZERLAND



AZ & PARTNER AG

Dorfchärn
CH-6243 Egolzwil
Switzerland

Tel. +41 41 982 02 02
Fax +41 41 982 03 03
www.azpartner.ch
info@azpartner.ch

AZ DENTAL GMBH

Siemensstrasse 5
DE-78467 Konstanz
Germany

Tel. +49 7531 50024

www.azdental.de
info@azdental.de

Distributed by: